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PEA 技术: olink.com/content/uploads/2021/09/olink-white-paper-pea-a-high-multiplex-immunoassay-technology-with-qpcr-or-ngs-readout-v1.0.pdf

Olink Explore 3072: olink.com/products-services/explore/

Olink Target 96: olink.com/products-services/target/#relative

Olink Target 48 Cytokine: olink.com/products-services/target/48-cytokine-panel/

Olink Flex: olink.com/products-services/flex/

Olink Focus: olink.com/products-services/custom-panels/

Olink NPX: olink.com/faq/what-is-npx/

定制型生物标志物panel的开发与验证: olink.com/content/uploads/2021/09/olink-development-and-validation-of-customized-pea-biomarker-panels-1083-v2.0.pdf

Olink Explore 和 Target panel的完整验证数据: olink.com/resources-support/document-download-center/

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1. Calderón-Celis F, et al. Standardization approaches in absolute quantitative proteomics with mass spectrometry. *Mass Spectrom Rev.* 2018 Nov;37(6):715-737.
2. WHO Expert Committee on Biological Standardization: report of the seventy-second and seventy-third meetings. Geneva: World Health Organization; 2021 (WHO Technical Report Series, No. 1030).
3. Neubert H, et al. Protein Biomarker Quantification by Immunoaffinity Liquid Chromatography-Tandem Mass Spectrometry: Current State and Future Vision. *Clin Chem.* 2020 Feb 1;66(2):282-301.
4. Palomares LA, Estrada-Mondaca S, Ramírez OT. Production of recombinant proteins: challenges and solutions. *Methods Mol Biol.* 2004;267:15-52.
5. Bratcher PE, Gaggar A. Factors influencing the measurement of plasma/serum surfactant protein D levels by ELISA. *PLoS One.* 2014 Nov 3;9(11):e111466.
6. Baskerville CL, et al. Comparison and validation of ELISA assays for plasma insulin-like growth factor-1 in the horse. *Open Vet J.* 2017;7(1):75-80.
7. WHO Expert Committee on Biological Standardization. Forty-third report. Geneva, World Health Organization, 2006 Annex 2 (WHO Technical Report Series, No. 932).
8. Enroth S, et al. High throughput proteomics identifies a high-accuracy 11 plasma protein biomarker signature for ovarian cancer. *Commun Biol.* 2019 Jun 20;2:221.
9. Brownridge P, Beynon RJ. The importance of the digest: proteolysis and absolute quantification in proteomics. *Methods.* 2011 Aug;54(4):351-60.
10. Lai X, Wang L, Witzmann FA. Issues and applications in label-free quantitative mass spectrometry. *Int J Proteomics.* 2013;2013:756039.
11. Olink®. Multiplex analysis of inflammatory proteins: A comparative study across multiple platforms. White paper. 2021 Sept.
12. Skau E, Wagner P, Leppert J. et al. Are the results from a multiplex proteomic assay and a conventional immunoassay for NT-proBNP and GDF-15 comparable?. *Clin Proteom.* 2023; 20, 5.
13. Zhao Y, Brasier AR. Qualification and Verification of Protein Biomarker Candidates. *Adv Exp Med Biol.* 2016;919:493-514.

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